

















DNA/RNA Sample Preparation Guide

	<p>Please write* the sample name clearly on the lid or side of the tube *Oil-based pen</p>		<p>Sample Volume : 25ul~100ul Sample Amount : Please prepare more than the amount that meet the QC criteria</p>
	<p>Frozen Shipping ※ Refrigerated delivery is possible if using the stable kit</p>		<p>Specific and accurate sample information is required. (Sample type, buffer, extraction methods & etc.,)</p>
	<p>1.5ml~2ml round tube or conical tube</p>		<p>If you use 96well plate, arrange the samples vertically</p>

DNA/RNA Extraction Sample Guide

 <p>Tissue</p>	<p>Over 50mg</p>	<p>DNA extraction 1) Add a RNAlater & frozen shipping</p> <p>RNA extraction 1) Add a RNAlater & frozen shipping 2) Add a Trizol & frozen shipping</p>	 <p>Cell</p>	<p>Over 2×10^6</p>	<p>DNA extraction 1) Cell pellet and add a media : delivery with frozen</p> <p>RNA extraction 1) Add a Trizol & frozen shipping 2) Add a media and delivery with frozen</p>
 <p>Blood</p>	<p>Over 200ul</p>	<p>DNA extraction 1) Using EDTA tube & refrigerated shipping</p> <p>RNA extraction 1) Use PAXgene or Tempus tube & refrigerated shipping 2) Using the EDTA tube & refrigerated shipping</p>	 <p>Plant</p>	<p>Over 50mg</p>	<p>DNA extraction 1) Add a RNAlater & frozen shipping</p> <p>RNA extraction 1) After freeze-drying, grind and frozen shipping</p>
 <p>Bacteria</p>	<p>Over 1×10^9</p>	<p>DNA/RNA extraction 1) Bacteria pellet : delivery with frozen 2) Culture with liquid media & delivery with refrigerated * The solid media(plate) is highly likely to be contaminated, so only refrigerated shipping is recommended.</p>	 <p>Serum/Plasma Exosome</p>	<p>Over 1ml</p>	<p>DNA extraction 1) Add a RNAlater & frozen shipping</p> <p>RNA extraction 1) Frozen shipping * Add 200ul elution solution to the sample from which Exosome is separated.</p>
	<p>If sample are not applicable to the guide, Please discuss with regional sales representative.</p>				

RNAlater / PAXgene Blood RNA tube Guide

Name of Kit		Maker	Relatives	
RNAprotect Tissue Reagents		QIAGEN (#76104)	<ol style="list-style-type: none"> 1. Add enough tissue to soak in the reagent. 2. Not recommend to use Cell, Blood, plasma and serum. 3. Not recommend for fatty tissues ¹⁾ 	<p>For RNA stabilization</p> <ul style="list-style-type: none"> ❖ 15~25°C, 1 weeks ❖ 2~8°C, 4 weeks ❖ -20 or -80°C, long-term storage
PAXgene Blood RNA tube		QIAGEN (#762165)	<ol style="list-style-type: none"> 1. After collecting blood, mix more than 10 times by inverting. 2. Store at RT for 2 to 72hrs and move to refrigerator (2-8°C) or Freezer(-20°C) ²⁾ 	<p>For RNA stabilization</p> <ul style="list-style-type: none"> ❖ 18~25°C, 3days ❖ 2~8°C, 5days ❖ -20 or -70°C, long-term storage
GenTegra RNAAssure™		GenTegra (#GTR50-LQ)	<ol style="list-style-type: none"> 1. Please use it as the final elution buffer for RNA extraction. 2. Recommend minimum elution volume is more than 20ul. ³⁾ 	<p>For RNA stabilization</p> <ul style="list-style-type: none"> ❖ 15~25°C, 3days ❖ 2~8°C, 2weeks ❖ -20 or -70°C, 1 year

1) <https://www.qiagen.com/us/products/discovery-and-translational-research/dna-rna-purification/rna-purification/total-rna/rnaproTECT-reagents-and-tubes/>

2) <https://www.qiagen.com/us/products/discovery-and-translational-research/sample-collection-stabilization/rna/paxgene-blood-rna-tubes/>

3) <https://gentegra.com/gentegra-rnassure/>